



Does your horse suffer from respiratory problems, allergies and airway irritations?

Would your horses improve their performance with clearer airways?

The Portable Nebulisation System for Horses

- **Simple & integrated**

Easy set up and operation no hoses, wires or awkward valves

- **Rapid**

10mls saline aerosolised in less than 12 minutes

- **Silent**

Noise Free

- **Adaptable**

For use with MDIs/inhalers

- **Savings**

On drug administration



fast, portable, equine nebulisation device

Worldwide distributors of the Flexineb™ equine nebuliser

**BreathEazy**®  
respiratory relief for cats, dogs and horses

For more information contact:

**BreathEazy**

154 Worcester Road, Malvern,  
Worcestershire WR14 1AA.  
United Kingdom

t: +44 (0) 845 680 8975

e: [enquiries@breatheazy.co.uk](mailto:enquiries@breatheazy.co.uk)

w: [www.breatheazy.co.uk](http://www.breatheazy.co.uk)



**BreathEazy**®  
respiratory relief for cats, dogs and horses



Flexineb™ is a fast, portable equine nebuliser device that is simple to use and noise free. It operates without the noise of loud compressors or the inconvenience of cumbersome hoses and pipes.

Flexineb utilises Nortev proprietary technology and is supremely portable allowing in excess of 3 hours of operation before the Lithium Ion battery needs recharging.

### Drugs:

The Nortev Flexineb system is designed to aerosolise a wide range of drugs including corticosteroids, bronchodilators, saline and antibiotics.

Medications must be prescribed by veterinarians; end users should follow the advice of veterinarians at all times.

### Cleaning & Disinfecting:

Cleaning is easy - simply wash with warm soapy water. Disinfecting is simple too, by using cold sterilisation methods.

### Volume of Medication:

A standard 10mls medication cup is accommodated. A 20mls extension is provided to facilitate volumes of up to 30mls. No medication remains in the cup after treatment.

### Comfort & Tolerance:


Flexineb fits comfortably onto the horse's nose as a result of the flexible material. An adjustable lever allows full control of the exit valve, enabling the caregiver to bypass the valve for easy breathing or to facilitate inhalation deep into the lungs.

### Visibility:

The transparent aerosol holding chamber allows the user to see the aerosol mist being generated by the Flexineb technology, in conjunction with the horse's breathing emptying the chamber and receiving the fine particle aerosol.

### Accessories:

System accessories include an AC/DC power supply to charge the internal battery, along with a user manual and a convenient storage bag.



1

Portable, easy to assemble mask



2

Battery operated controller



3

Nebuliser