

The Portable Nebulisation System for dogs, cats and other small animals

#### Silent

Completely silent when in use

#### Simple to Use

Minimal number of parts and a one touch on/off button

### Portable

Battery operated so can be used anywhere

#### Fast

Faster treatment times compared to compressor nebulisers

### **Cost Effective**

Cost savings on drugs compared to other therapies such as MDIs (Meter Dose Inhalers)



For more information contact: BreathEazy Ltd 154 Worcester Road, Malvern, Worcestershire WR14 1AA United Kingdom

t: +44 330 053 5056 f: +44 330 043 5056 e: enquiries@breatheazy.co.uk w: www.breatheazy.co.uk



Fast, Portable, Versatile Small Animal Nebulisation Device



# 

#### What is Nebulisation?

Nebulisation is a treatment technique in which tiny spherical droplets of aerosolised medication are generated and inhaled by the patient animal.

#### Why use a Nebuliser?

A nebuliser device enables the small animal owner to easily administer care to their pets' airways and lungs. The device incorporates a mask that fits on the nose and a technological means of generating a fine mist.

### Why Flexineb® c1

- Flexineb<sub>C1</sub> is a fully integrated nebulisation system designed to treat a number of respiratory diseases in small animals including cats, dogs, rabbits, rodents, reptiles and birds.
- Flexineb<sub>C1</sub> incorporates the same vibrating mesh nebulisation technology used in Nortev's award winning Flexineb<sup>™</sup> equine nebuliser.

- Flexineb<sub>C1</sub> can be used to treat a number of conditions such as feline asthma and allergic rhinitis in cats; and chronic bronchitis and laryngeal paralysis in dogs. Flexineb<sub>C1</sub> can also be used to treat a number of respiratory infections associated with small animals.
- You can use a wide range of drugs with Flexineb<sub>C1</sub> including antibiotics, bronchodilators, steroids and mucolytics. You can also use non-drug based solutions such as saline.
- The Flexineb<sub>C1</sub> has 2 types of interface to treat your pet. The first uses soft silicone masks as the interface between the nebuliser and patient. This makes the treatment portable. These silicone masks come in 3 sizes, (small, medium and large) to suit all patient types.
- The second interface is a small animal aerosol treatment chamber. This can be used for cats and small dogs that are unwilling to be treated using the silicon mask interface or in an emergency situation. The treatment chamber can also be used to treat a number of other small animals including rabbits, guinea pigs, birds and tortoises.



## Portable, easy to assemble system

# Battery operated controller

# Medication cup and handheld spacer

